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

*Precise temperature control
Fast heating rates
Little space requirement
Horizontal or vertical operation
Arbitrary working gases*

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The LORA type furnaces produced by HTM Reetz GmbH are very versatile models.


*This type of furnace is capable of working up to 1800 °C in inert, reducing or oxidizing atmosphere
fast heating rates possible
with only little space requirement
can be operated in any position*

For specialist applications the LORA furnace may be designed as gradient or multiple-zone tube furnace.

On demand, its external case may be equipped with a water cooling capability.

For operation in vacuum and for other specialist applications, we offer a furnace version with insertion tube.

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General Design

The essential part of the furnace is the work tube made of sintered aluminum oxide. The molybdenum heating element wire is directly wound on the work tube. The workspace of the furnace is formed by the interior of the tube.

The space outside the tube is the heating element area which is hermetically encased into a stainless steel housing.

The standard dimensions of the work tube are: interior diameter 32 mm, length 600 mm. If technically feasible, other dimensions are available.

The Workspace

ISO-KF standard vacuum flanges close up the ends of the workspace. Swagelok® connectors for gas inlet and outlet are integrated into the flanges which allow flushing the workspace with any gases. Within limits, work in vacuum is possible.

The temperature is measured by a type B thermocouple covered by a protection tube. It is placed in the workspace and fed outside the furnace through a flange.

The Heating Element Area


The stainless steel housing encasing the heating conductor space contains the thermal insulation as well as a getter that protects the heater from oxidation. During operation, the heating element area of the tube furnace must be flushed with a slight inert gas flow.

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
Temperaturprofile / Temperature profiles


Standard-Rohröfen
Standard tube furnace

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Separate temperature control units

The temperature control units for all LORA tube furnaces offered by HTM are designed as separate units in 19" cases. The indicating and operating displays of the controller, the load switch and main switch, the ammeter are placed on the front panel. The mains and load connections as well as the connections for thermocouple and serial interface are on the back panel.

PID programming EUROTHERM or OMRON controllers are used in the temperature control unit. Power control is effected by phase-controllable thyristors.

The control unit is equipped with a relay module for automatic cut-off in case of excess temperature. Optionally it can be equipped with RS 232, RS 422 or RS 485 interfaces.

Based on an additional controller and separate thermocouple an independent temperature control can be integrated in the control unit.

Software "HTMCtrl+" enables optimal operation

Combined with our software "HTMCtrl+" the operation of the controllers can be simplified and a simple programming controller can be upgraded to a multi-programming controller with much more convenience. It allows simultaneous communication with max. 7 controllers.

"HTMCtrl+" software enables:

- configuring and programming of controllers*
- storing of data and programs*
- graphic and numeric visualization of data*
- graphic display of the total period as well as partial periods of an experiment*
- storing of measuring data in EXCEL format CSV*
- control of max. 4 MFC in correlation with the heating program*

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Modular design

The temperature profile of the furnace can be specifically influenced by varying the density of the heating wire winding. As a result HTM can offer for example furnaces with a long zone of uniform temperature gradient furnaces

Due to the technical configuration and the variable working position of the LORA furnace the application possibilities can be extended by using additional modules as for example

*insertion tube, e.g. for vacuum operation
vacuum pump
thermal balance
window flange
mass flow controller
charging device
oxygen sensor*

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HTM Reetz GmbH, USA